COMPARISON OF RARE GERMLINE VARIANTS ASSOCIATED WITH CANCER PREDISPOSITION SYNDROMES IN TAIWAN AND USA POPULATIONS

Zulfiqar Gulzar¹, Julian Lucas¹, Alex Atkins¹, Stephen Su¹, Oscar Segurado², Rui Mei¹*

1. CellMax Life, 1271 Oakmead Parkway, Sunnyvale, CA 94085

2. MedicAffairs Consulting LLC, 6762 Hampton Dr, San Jose, CA 95120

* Corresponding author: rui@cellmaxlife.com

BACKGROUND AND PURPOSE

Germline mutations that increase the susceptibility of an individual to certain cancers have been identified in numerous genes, and individuals can be tested for mutations in these genes to understand their inherited lifetime risk for developing cancer ¹⁻⁵. Traditional approaches, using Sanger sequencing or qPCR, focus on small groups of genes and therefore present a difficulty to screen for mutations in numerous genes from multiple individuals synchronically. As the landscape of genes that are important to hereditary cancer risk grows, it is becoming increasingly important to develop methods to interrogate multiple genes that may carry mutations. With rapid advances in next-generation sequencing (NGS) technology and bioinformatics, testing of multiple genes from multiple individuals synchonically with panelbased tests is becoming increasingly common.

The CellMax DNA Genetic Cancer Risk Test was designed to allow the analysis of genes and genes greatly affects cancer risk. genetic variants associated with predisposition to cancer. The test targets germline variants detectable in white blood cell DNA. The 98 gene panel [Fig.1] was selected based on systematic, expert review of the scientific literature to identify genes that predispose to cancer, and focuses on genes for which there is already strong evidence of such predisposition. The focus of this study was to evaluate the analytical and clinical performance of the CellMax DNA Genetic Cancer Risk Test, and to observe the prevalence of highly-penetrant, rare, pathogenic germline variants strongly linked to predisposition of cancers in patient populations from Taiwan and USA.

METHODS AND STUDY DESIGN

The CellMax DNA Genetic Cancer Risk Test covers all coding exons, nearby flanking regions (+/- 20bp), and known splice-sites in the 98 gene panel. The test detects single nucleotide substitutions, insertions and deletions in genomic DNA with a mean coverage of >100X for each sample.

As part of this study, we validated the 98 gene panel prior to sequencing the patient samples from the study cohorts. The assay validation design completely took the advantage of the availability of gold standard reference samples. The lab performed end-to-end sample processing and data analysis with 19 unique gold standard reference samples run over 38 different library preprations throughout the analytical validation. After processing from DNA input through variant calling [Fig.2], the data produced is compared to the high-quality reference calls to assess the assay's analytical performance. A concordance study to compare genotyping performance between blood and saliva-derived DNA was also performed.

We collected **2,259 patient samples** from two clinical centers for testing on the CellMax DNA Genetic ready libraries, enriched for target regions, sequenced, and variants are called Cancer Risk Test 98 gene panel - 1885 patients from CellMax's CAP laboratory in Taiwan, and 374 patients from CellMax's CAP/CLIA laboratory in California, many of whom had been referred for hereditary cancer risk counseling or assessment. Blood or saliva specimens were accepted for sample processing. Information on demographics, personal, and family history of cancer was collected at time of sample collection. The mean age of the Taiwan cohort was 43, with 58% males and 42% females, while the US cohort had a mean age of 58 with 40% males and 60% females.

Sequence variant classification as pathogenic, likely pathogenic, variant of uncertain significance (VUS), likely benign and benign was performed according to the American College of Medical Genetics and Genomics (ACMG) sequence variant interpretation guidelines ⁶. All classifications were evaluated by a board-certified pathologist and only pathogenic variants were clinically reported.

REFERENCES:

Abstract 2241 / 23 AACR Annual Meeting 2018 Chicago, IL

Expert-C	Expert-Curated, Comprehensive 98-Gene Panel							
AIP	ALK	APC	ATM	BAP1	BARD1	BLM	BMPR1A	BRCA1
BRCA2	BRIP1	BUB1B	CDC73	CDH1	CDK4	CDKN1C	CDKN2A	CEBPA
CEP57	CHEK2	CYLD	DDB2	DICER1	DIS3L2	EGFR	EPCAM	ERCC2
ERCC3	ERCC4	ERCC5	EXT1	EXT2	EZH2	FANCA	FANCB	FANCC
FANCD2	FANCE	FANCF	FANCG	FANCI	FANCL	FANCM	FH	FLCN
GATA2	GPC3	HNF1A	HOXB13	HRAS	KIT	MAX	MEN1	MET
MLH1	MSH2	MSH6	MUTYH	NBN	NF1	NF2	NSD1	PALB2
РНОХ2В	PMS1	PMS2	PPM1D	PRF1	PRKAR1A	PTCH1	PTEN	RAD51C
RAD51D	RB1	RECQL4	RET	RHBDF2	RUNX1	SBDS	SDHAF2	SDHB
SDHC	SDHD	SLX4	SMAD4	SMARCA4	SMARCB1	STK11	SUFU	TMEM127
TP53	TSC1	TSC2	VHL	WT1	WRN	XPA	ХРС	

[Fig.1] The 98 genes included in the CellMax Life DNA Genetic Cancer Risk Test. Genes associated with hereditary cancer dromes and cancer risk were included when there was sufficient evidence in the literature that mutations in thos



[Fig.2] In the CellMax DNA Genetic Cancer Risk Test workflow, genomic DNA is converted to sequence-



RESULTS

The analytical sensitivity, specificity and accuracy of the assay panel were 99.95%, 100% and 99.99%, respectively. Similarly, the reportable variants for 30 DNA samples containing clinically important SNVs in BRCA1/2, RET, APC, MEN1 and PTEN were called with no false positives. Paired blood (Reference lab) and saliva (CellMax CLIA lab) samples were run to validate the sample collection protocol. The concordance for the paired blood/saliva sample analysis was found to be **99.99%**.

Summary of Results

• Germline pathogenic variants were detected in 47/1885 (2.5%) and 25/374 (6.7%) of individuals in the Taiwan and US cohort respectively [Fig.4].

• Pathogenic BRCA1 and BRCA2 mutations comprised 8.5% and 17% of total positive cases in Taiwan, and **4% and 0%** of total positive cases in the US.

• Outside of BRCA1 and BRCA2, additional cases of rare pathogenic mutations in Taiwan 35/47 (75%) and the US 24/25 (96%) associated with increased cancer risk were identified with our panel [Fig.5], which would not have been identified if only BRCA1/2 testing was used. These genes include MUTYH, MLH1, RET, ATM, APC, CHEK2, ERCC2, PALB2, PTCH1, WRN, BARD1, BRIP1, PMS2, PTEN, RAD51D, SDHD, RUNX1, BLM, HOXB13, and MSH6 [Fig.6].

• 76% of non-BRCA1/2 findings in the US and 40.5% of non-BRCA1/2 findings in Taiwan are in genes which have well-defined national guidelines outlining care and would warrant consideration of a change in patient care.

• The pathogenic variants within non-BRCA1/2 genes are associated with breast/ovarian cancer genes (ATM, CHEK2, PALB2), MUTYH-associated Polyposis syndrome gene (MUTYH), and Lynch syndrome genes (MLH1, MSH2, MSH6, PMS2).



were more common in the US cohort

CONCLUSIONS

We have demonstrated in 2 clinical cohorts the analytical and clinical validity of a 98 gene panel for the detection of rare pathogenic germline variants associated with cancer predisposition. High sensitivity and specificity of this panel was demonstrated irrespective of the use of blood or saliva as samples. This study provides useful data of the prevalence for many cancer predisposition genes in two diverse ethnic groups. The frequencies of pathogenic mutations within these cohorts revealed a disparity between racial/ethnic groups with similar risks. Multi-gene panel testing using next-generation sequencing provides more comprehensive, cost effective, and clinically actionable assessment compared to traditional single gene hotspot testing for hereditary cancer risk.

CellMax Life

AACER American Association for Cancer Research

ANALYTICAL VALIDATION SUMMARY

Analytical Validation	Variant Present	Variant Not Present	
Variant Detected	1841 true positives	0 false positives	
Variant Not Detected	1 false negative	3,898,611 expected reference calls	
Results	99.95% Sensitivity	100% Specificity	

[Fig.3] Three sets of gold standard reference samples and data run over several different library preparations were used to analytically validate the end-to-end sample process and data analysis.

PERCENTAGE OF PATHOGENIC VARIANTS IN TWO DIVERSE CLINICAL POPULATIONS





[Fig.4] The study cohorts consisted of 374 patients from the US and 1885 patients from Taiwan. These patients were referred for hereditary cancer risk counseling or assessment. Prevalence of pathogenic mutations was observed to be higher in the US cohort.



[Fig.5] BRCA1/2 testing alone would only have identified 12 individuals from the Taiwan cohort and 1 individual in the US cohort gene panel testing identified other important pathogenic mutations in 35 additional patients from Taiwan and 24 onal patients from the US.